


# Fermentation Lab Presentation

**ABE 304**

# Peer Evaluations

- 4 peer assessments throughout the semester
  - **Open February 6 - February 12 (42% response rate)**
  - Open February 24 – February 28
  - Open March 23 – March 28
  - Open April 27 – May 2
- Administered through CATME

# Fermentation



# Fermentation Products


- How many fermentation products can you think of?

# Fermentation Products


- Ethanol
  - Beer
  - Wine
  - Liquor
- Yogurt
- Kombucha
- Kefir
- Sauerkraut
- Kimchi
- Pickles
- Cheese
- Salami
- Prosciutto
- Sourdough bread
- Pharmaceutical products
  - Penicillin
  - Insulin
- Citric Acid
- Acetic Acid

# Xanthan Gum

- Microbial Polysaccharide
- Produced from *Xanthamonas campestris*
- Significant rheological properties
- Cannot be synthesized



# Industrial Production



# Industrial Fermentation

- Temperature uniformity
- Reactor size
- Internal pressure
- pH
- Oxygen distribution
- Nutrient availability
- Level of agitation




# In the Lab


- Study effect of temperature on cell growth and xanthan production

25°C	30°C
Groups 1 & 2	Groups 3 & 4

# Growth Kinetics



# Growth Kinetics



$$r_x = \mu x \qquad r_x = \frac{dx}{dt}$$


$$\frac{dx}{dt} = \mu x$$

$$\ln x = \ln x_0 + \mu t$$

$$x = x_0 e^{\mu t}$$

# Lab – Week 1

1. Inoculate fermentation medium with *Xanthamonas campestris*
2. Take samples every 45 minutes for a total of 4 samples
3. Take samples from previously inoculated media to capture all phases of growth (12 samples total)



# Cultures

- Each group will have three cultures
  - First—Two days prior to lab (TA) → 4 samples
  - Second—One day prior to lab (TA) → 4 samples
  - Third—Inoculate *Xanthamonas campestris* in class (Group) → 4 samples

5mL inoculum → 250 mL flask (with 100mL media)




Incubate your flask in the shaker at 25°C or 30°C (~150rpm)

# Yield

- Measure of how much of the substrate is converted into product
- Allows us to better account for
  - Reaction efficiency
  - Side Reactions


# Calibration

Each spectrophotometer must be calibrated to be able to report concentration




# Dilution

Some samples will need to be diluted to obtain accurate results



= 1 mg/mL

Add 9 mL of water



= 0.1 mg/mL

How  
diluted is  
the  
sample?



# Lab – Week 1

## Cell Growth



☐ Take sample

☐ Cell density (Spectrophotometer)

☐ Centrifuge

☐ Remove supernatant

☐ Measure with Spectrophotometer (at 600 nm)

Blank: Media before inoculation

## Supernatant



☐ Save supernatant in freezer for week 2

# Sample

- 1 mL broth before inoculation (reference blank)
- 1 mL sample for Cell Density
- 1 mL sample x 2 = 2 supernatant sample for Glucose & Xanthan concentrations
  - Supernatant will be frozen for use during week 2

Total volume of each sample = 3mL

# Lab – Week 2

## Glucose Concentration



☐ Use frozen supernatant and dilute with water

☐ Add glucose oxidase (GOPOD kit)

☐ Measure with spectrophotometer (at 510 nm)  
Blank: GOPOD kit solution

## Xanthan Concentration



☐ Xanthan calibration curve ( Xanthan stock solution)

☐ Use frozen supernatant and dilute with water

☐ Sulfuric Acid digestion

☐ Measure with spectrophotometer (at 325 nm)  
Blank: Water