

Bioprocess Engineering Lab: Separations in Bioprocessing

ABE 304

Why are Separations important?

- Biological components are rarely found alone

Why are Separations important?

- Biological components are rarely found alone
 - Biomass components
 - Pharmaceuticals
 - Proteins
 - Bacteria
 - Yeast

BioSeparations... It's a process



1. Separate Large solids


2. Remove Smaller
Particles

3. Purify
Desired
Component(s)

BioSeparations... It's a process

- Solid/Liquid
 - Sedimentation
 - Centrifugation
 - Filtration
- Recovery
 - Membrane
 - Precipitation
 - Extraction
- Purification
 - Crystallization
 - Chromatography

Unit Operation	Level at which separation occurs
Centrifugation	Macroscopic
Filtration	Macroscopic
Microfiltration	Macroscopic
Ultrafiltration	Molecular
Reverse Osmosis	Molecular
Crystallization	Macroscopic
Extraction	Molecular
Adsorption Desorption	Molecular Molecular
Chromatography	Molecular
Affinity Chromatography	Molecular




Sedimentation

- Gravity driven separation
 - Large particle size
 - Difference in density of solid from the liquid
- Common in wastewater applications




Sedimentation

- Gravity driven separation
 - Large particle size
 - Difference in density of solid from the liquid
- Common in wastewater applications
- Also used to separate Brewers yeast
 - CO₂ created during fermentation can complicate the settling
 - Cells settle out once fermentation is finished



Sedimentation


- Can manipulate gravity to decrease separation time



Centrifugation

- Enhanced sedimentation
 - Multiply the force of gravity by rotating the container
 - Centrifugal force function of diameter and rotational speed
 - Results in two streams
 - Clarified liquid (supernatant)
 - Wet solids

Centrifugation




Centrifugation

- Batch process



- Continuous Process


Disk centrifuge



https://www.youtube.com/watch?v=zSL_-DcsjR4


Filtration

- Used when solids are large enough to be stopped by a filter cake
 - Either pumped



Filtration


- Used when solids are large enough to be stopped by a filter cake
 - Either pumped or **pulled** through



Membrane Separations

- Membranes enable filtration of
 - Colloids
 - Cells
 - Molecules
- By means of
 - Microfiltration – remove microorganisms to a 0.45 or 0.22 μm cut off
 - Ultrafiltration – cut off of 5 to 5000nm
 - Reverse Osmosis – semipermeable with a pore size of 0.5 to 5 nm


Membranes



Note: 1 Micron (1×10^{-6} Meters) = 4×10^{-5} Inches (0.00004 Inches)
1 Angstrom Unit = 10^{-10} Meters = 10^{-4} Micrometers (Microns)

© Copyright 1996, 1993, 1990, 1984 Osmonics, Inc., Minnetonka, Minnesota USA




Membranes



Microfiltration


- Organic or inorganic membranes
 - Pores 0.1 to 10 microns
- Membrane structure
 - Flat, pleated sheet
 - Hollow fiber

Membrane Filtration




Microfiltration

- Membrane structure
 - Hollow fiber
 - Individual fibers are 1 to 2 mm in diameter
 - Bundled together for more surface area



Precipitation




The other kind of precipitation

- Salting-out
 - Hydrophobic
 - Decreases the solubility of the solute in solution
- Salting-in
 - Electrostatic
 - Increase solubility with increasing salt
- Example: [Proteins in milk](#)

Extraction

- Uses solubility differences of biological molecules between
 - Two different liquid phases
 - A solid and a liquid



Soxhlet Extractor

Crystallization

- Nucleation process that occurs in a concentrated solution
- Nucleation occurs when
 - The solute concentration exceeds the saturation AND
 - Solute molecules are present in large enough clusters to serve as nucleation sites



Supersaturation

- For crystals to form, the solution must be super saturated
- How to get a super-saturated solution?
 - Cool a solution
 - Works if the solubility of the salt in the solution decreases with temperature
 - Evaporate a portion of the solution
 - Add a third component
 - Salting-out


Chromatography

- Partitioning effects based on
 - Differential binding
 - Molecular sieving
- Partitioning of the target molecule relative to other molecules in the solution
- Stationary phase binding
 - Ion exchange
 - Reversed and normal phase
 - Hydrophobic interaction
 - affinity

Purification Types

- Stationary phase binding
 - Ion exchange
 - Reversed and normal phase
 - Hydrophobic interaction
 - Affinity
- Size Exclusion
 - Gel permeation
 - Partitions with a sieving effect
 - Small molecules penetrate the stationary phase
 - Larger molecules are excluded and elute first


Chromatography Nomenclature



Frits are sintered metal or polymeric filters with particle cut-off from 2 to 20 microns




Column is packed with spherical particles (5 – 300 micron)



Chromatography Nomenclature

Analytical scale are operated either incline, vertical or horizontal

Process scale columns are usually vertical



Less strongly retained solute elutes before the more strongly retained solute.

BioSeparations... It's a process



1. Separate Large solids

2. Remove Smaller
Particles

3. Purify
Desired
Component(s)

Xanthan Gum purification and concentration

BIOSEPARATIONS LAB

Read through Chapter 11 in *Bioprocess Engineering Principles* for additional help

Types of Membrane Filtration





Figure 11.23 Doran

Process for Lab



Phenomena to Explore

- Retention coefficient $R = \frac{C_R - C_P}{C_R}$
- Permeate Flux $J = \frac{F_p}{A}$
- Mass Transfer Coefficient $J = k \ln \frac{C_F}{C_B}$

Control Regime




Figure 11.2 - Doran

Use your data to determine if filtration is pressure-controlled or mass-transfer controlled.




Figure 4.11 - Cheryan, 1998
Ultrafiltration and Microfiltration Handbook

Similar to Figure 11.30 - Doran

Basic Lab Procedures

- Calibrate the spectrophotometer
- Calibrate the peristaltic pump
- Test 4 flow rates at 5 transmembrane pressures (20 pairs of conditions)
- At each setting, collect samples from the permeate and retentate (40 samples total)
- Determine xanthan concentrations using spectrophotometer